

## **Total Ash in Forages**

### **1. Application**

This procedure is applicable for the determination of ash in all types of dried, ground forages and feeds.

### **2. Summary of Methods**

A dried, ground sample is ignited in a furnace at 500°C to oxidize all organic matter. Ash is determined by weighing the resulting inorganic residue.

### **3. Safety**

Basic precautions regarding mechanical equipment, electric motors, and glassware must be followed. All electrical equipment is properly grounded and installed and maintained by qualified electricians.

### **4. Interferences**

This procedure is not applicable for ash determination in liquid feeds or feeds high in sugar content.

### **5. Sample Collection, Preservation, and Handling**

All samples are dried at 55°C in a cabinet-type forced air dryer for 12-16 hours. After drying the sample is ground to pass through a 1 mm forage mill. A subsample is then dried at 105°C for 3 hours to determine laboratory DM content.

### **6. Apparatus and Materials**

- 6.1 30 ml porcelain crucibles, low wide form, numbered with furnace-proof ink
- 6.2 Muffle furnace with pyrometric controller
- 6.3 Analytical balance, sensitive to 0.1 mg
- 6.4 Desiccator, with vented lid
- 6.5 Drying oven

### **7. Reagents**

None

## 8. Methods

- 8.1 Remove crucibles, which have been dried for at least 2 hours at 105°C from oven, and move to desiccator. Cool and record weight of crucibles to the nearest 0.1 mg ( $W_1$ ).
- 8.2 Weigh 1.5 to 2.0 g of sample into the crucible, recording weight of crucible and sample to the nearest 0.1 mg ( $W_2$ ).
- 8.3 Ash in furnace at 500°C for 2 hours after the furnace reaches temperature.
- 8.4 Allow crucibles to cool in furnace to less than 200°C and place crucibles in desiccator with vented top. Cool and weigh crucible and ash to the nearest 0.1 mg ( $W_3$ ).

## 9. Calculations

- 9.1 % Ash (DM basis) =  $[(W_3 - W_1) * 100] / [(W_2 - W_1) * \text{lab DM} / 100]$   
 $W_1$  = tare weight of crucible in grams  
 $W_2$  = weight of crucible and sample in grams  
 $W_3$  = weight of crucible and ash in grams

## 10. Quality Control

An in-house standard is run to gauge run acceptability.

## 11. Reporting

Results are reported as % ash on a dry matter basis.

## 12. References

- 12.1 Ash of Animal Feed. (942.05) Official Methods of Analysis. 1990. Association of Official Analytical Chemists, 15<sup>th</sup> Edition.