

Organic Matter

Weight Loss-on-Ignition (LOI 360°)

1. Application

This procedure is used for the routine estimation of soil organic matter by the loss of weight in a sample heated at a temperature high enough to burn organic matter but not so high as to decompose carbonates.

2. Summary of Methods

A sample of soil is dried at 105° C to remove moisture. The sample is weighed, heated at 360° C for 2 hours and weighed again after the temperature drops below 150° C.

3. Safety

Care should be exercised in handling hot samples. Be sure to cool the oven to 150° C before removing the samples from the oven. Use a good pair of tongs and grasp the sample firmly.

4. Interferences

Any material that losses moisture below 360° C is a potential source of error. Therefore, soil moisture must be removed before the base weight of the sample is taken. Also, ignited samples must not be allowed to re-absorb moisture from the air before they are weighed.

Gypsum loses water of hydration gradually. Soils containing gypsum should be heated initially at 150° C instead of 105° C. Some hydrated clays may also lose water below 360° C.

It is important that the results of this method be calibrated against organic carbon, preferably using a carbon analyzer, on soils from the area for which the test will be used.

5. Apparatus and Materials

- 5.1 Oven, or muffle furnace capable of being heated to 400° C and controlled to within $\pm 10^\circ$ C.
- 5.2 Beakers, 20 ml
- 5.3 Crucible rack, stainless steel
- 5.4 Balance accurate to ± 0.001 g in a draft free, low humidity environment
- 5.5 Soil scoop calibrated to hold 5 g of light-colored silt loam soil
- 5.6 Drying oven, 105° C

6. Reagents

An advantage of this method is that no reagents are required.

7. Methods

- 7.1 Place a 5 g scoop of soil into a tared 20-ml beaker
- 7.2 Dry for 2 hours or longer at 105° C
- 7.3 Record weight to ± 0.001 g
- 7.4 Bring oven to 360° C. Samples must then remain at 360° C for two hours.
- 7.5 Cool to $< 150^\circ$ C
- 7.6 Weigh to ± 0.001 g, in a draft-free environment

8. Calculations

- 8.1 Calculate percent weight loss-on-ignition (LOI)

$$\text{LOI} = \frac{(\text{wt. at } 105^\circ\text{C}) - (\text{wt. at } 360^\circ\text{C})}{\text{Wt. at } 105^\circ\text{C}} \times 100$$

- 8.2 Estimate % organic matter. Organic matter is estimated from LOI using regression analysis. Select soils covering the range in organic matter expected in the area serviced by the lab. Determine % organic matter using a carbon analyzer or by the Walkley-Black procedure for organic carbon. Regress OM on LOI.

9. Quality Control

- 9.1 At least one standard soil of known LOI value should be run with each batch of samples. If the result is not within the known standard deviation, corrective action is required.
- 9.2 All beakers should be re-tared monthly. Two beakers from each batch of 50 should be re-tared weekly. If the results are not within ± 0.002 g of the previous tared weight; re-tare all beakers in the batch.

10. Reporting

Data are reported as % LOI or as estimated % O.M.

11. References

- 11.1 Combs, S.M., and Nathan, M.V. 1998. Soil organic matter. Pp. 57-58. *In* J.R. Brown (Ed.), Recommended Chemical Soil Test Procedure for the North Central Region. NCR Publ. N0. 221 (revised). Missouri Agr. Exp. Sta. SB 1001. Columbia, MO.
- 11.2 Schulte, E.E., and Hopkins, B.G. 1996. Estimation of soil organic matter by weight

loss-on-ignition. Pp.21-31. *In* F.R. Magdoff, M.A. Tabatabai, and E.A. Hanlon, Jr. (eds.), *Soil Organic Matter: Analysis and Interpretation*. Soil Sci. Soc. Am., Madison, WI.